For use with spinit® CRP test and spinit® instrument

**PRODUCT DESCRIPTION**

**Intended Use**  
spinit® CRP Control kit contains liquid preparations of delipidized human serum at two concentration levels (Control L1 and Control L2) of C-Reactive Protein (CRP). The controls should be used to confirm that your spinit® instrument system is working properly and provides reliable results. Only when controls are used routinely, and the values obtained are within acceptable ranges can accurate results for patient samples be assured.

**Composition**  
• 3 x 250 µL spinit® CRP Control L1  
• 3 x 250 µL spinit® CRP Control L2  
• Instructions for Use

**WARNINGS AND PRECAUTIONS**

- For in vitro diagnostics use.
- Do not use the spinit® CRP controls after the expiry date or if the controls have not been stored in accordance to recommendations.
- Proper handling and disposal methods should be followed in accordance with the applicable country laws;
- All individual human blood donor units used to manufacture this lot were tested by FDA approved methods and found to be nonreactive for HbsAG, HIV-1, HIV-2, and HCV.
- Handle and dispose of the spinit® CRP controls as potentially bio hazardous materials. Use gloves.

**STORAGE AND HANDLING**

- The expiry date of the kit only applies for unopened product stored refrigerated (2-8°C).
- Avoid exposure to direct sunlight.
- Do not freeze.

**Unopened control vials**

- Unopened control vials are stable until expiry date indicated on the vial label when stored refrigerated (2-8°C).

**Opened control vials**

- Opened control vials are stable for 30 days when stored refrigerated (2-8°C). It is recommended to note the date of opening on the vial label.
- Always store the control vials in an upright position at 2-8°C when not in use.

**FREQUENCY OF CONTROL TESTING**

Testing of quality controls should be performed as per the established laboratory program. biosurfit also recommends testing controls:

- When installing the spinit® instrument for the first time;
- When testing each new lot of spinit® CRP test discs;
- Any time an unexpected value test result is obtained;
- After installing a software update;
- When training new operators on how to use spinit® instrument;
- In compliance with national or local regulations.

**TESTING PROCEDURE**

- Confirm the instrument is set to QC mode when running controls.
- Allow the control material to reach room temperature (15-25°C) before use.
- Mix the control material well before collecting a sample. Do not invert the tube for mixing.
- The sample can be extracted directly from the vial using a spinit® sample strip.

 Consult the spinit® CRP instructions for Use.

**Analysing a control**

 Consult the spinit® instrument User manual or the spinit® CRP Quick Guide for a step-by- step procedure.

**Verifying the control results**

The measured values should be within the acceptable range stated for the spinit® CRP control L1 and L2 (consult label on each kit box).

If the result obtained is outside the acceptable limits, verify that:

- Patient samples are not analysed until control results are within acceptable limits;
- The control vial has not passed its expiry date;
- The control vial has not been in use for more than 30 days;
- The control vial and the spinit® CRP test cartridges have been stored according to recommendations;
- There is no evidence of bacterial or fungal contamination of the control vial;
- The controls were tested in QC mode.

Correct the procedural error and retest the control material.

If no procedural errors are detected:

- Retest the control material using a new control vial;
- Examine the laboratory's quality control record to investigate the frequency of control failures;
- Patient results must be declared invalid when controls do not perform as expected;
- Contact your local supplier for advice before analysing patient samples.

**Target values**  
The spinit® CRP test used for establishing the specification for spinit® CRP Controls is traceable to the ERM™-DA472/IFCC2 reference. Target values and acceptable ranges for the spinit® CRP controls are stated on the kit box.

**SYMBOLS**

				
Manufacturer's name and address	Consult instructions for use	Lot number	Date to use by	In vitro diagnostic medical device

		
CE marking of conformity to directive 98/79/CE on in vitro diagnostic medical devices	Storage temperature between 2-8°C	Biological risk

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